

PC Catalog



**Panda 🐾
Scanner**

Changes in this version of the PC catalog

[illegible]

The recommended configuration for the Panda Scanner series is as follows






Recommended computer configuration	
CPU	Intel i7-12700H; Intel i7-13650H; Intel i7-13700H; Intel i7-14650HX
RAM	16G at least, 32G for better performance
Hard disk	SSD 512G or above
GPU	NVIDIA GeForce RTX 3060、4060, with at least 6G memory
Operating system	Windows10/11 64bit
NVIDIA Driver	Version 536 or above

- 1 The computer needs an independent graphics card. Only NVIDIA graphics cards are supported.
- 2 AMD processor is not supported.
- 3 MacOS is not supported.
- 4 Not recommended brand: HP; ASUS

Minimum computer configuration	
CPU	Intel i7-11800H
RAM	16G
Hard disk	512GB
GPU	RTX 2060
Operating system	Windows10/11 64bit



Here are some computers we recommend

Item		PANDA smart/PANDA free	PANDA P3/P4
Dell Alienware X16 CPU i9-13900HX GPU RTX 4060 RAM 32GB RAM DISK 512GB SSD		✓	✓
Dell G15-5530 CPU i7-13650HX GPU RTX 4060 RAM 16GB DISK 512GB SSD		✓	✓
Lenovo IdeaPad 3i CPU i7-12650H GPU RTX 3060 RAM 16GB DISK 512GB SSD		✓	✓
MSI GL76 CPU i7-12650H GPU RTX 3060 RAM 16GB DISK 512GB SSD		✓	✓
Dell G15 5520 CPU i7-12700H GPU RTX 3060 RAM 16GB DISK 512GB SSD		✓	✓

For more details, please visit:

<https://www.amazon.com/>


<https://www.lenovo.com/>

<https://www.dell.com/>

How to Check Your Computer Information

You can check your computer information using the following method:

Open **Settings** > **System** > **About**.



Control Panel > **System** > **Device manager** > **Display adapters**

