PC Catalog



www.panda-scanner.com info@panda-scanner.com +86 512 6731 7873


Panda
Scanner

Changes in this version of the PC catalog

Date	Details
2023/8/15	Update Configuration

The recommended configuration for the Panda Scanner series is as follows

Recommended computer configuration	
CPU	Intel i7-12700H/Intel i7-11800H or above
RAM	16G/32G
Hard disk	SSD 512G or above
GPU	RTX2060/RTX3060 or above
Operating system	Windows10/11 64bit

1. The computer must be equipped with an independent graphics card, and it is necessary to ensure that the memory of the graphics card is more than 6GB, and our software does not support AMD graphics cards at present.

2. The CPU for Intel is recommended, but the AMD is not support.

3. We recommend using a laptop designed for gaming or using a workstation, but not a light or office laptop. Computers not powered by Windows such as Mac is not support too.

4. At the same time, Panda Smart has higher CPU requirements and we recommend using I7 12th generation or more advanced CPUs to achieve a best performance.



Here are some computers we recommend

Item		PANDA smart	PANDA P2/P3
Alienware m16 Gaming Laptop CPU: i9-13900HX RAM: 16 GB, 2 x 8 GB, Disk: 1 TB SSD GPU: RTX 4080 12 GB		✓	✓
Alienware m16 Gaming Laptop CPU: i7-13700HX RAM: 16 GB, 2 x 8 GB Disk: 1 TB SSD GPU: RTX 4070 8 GB		✓	✓
G16 Gaming Laptop CPU: i7-12700H RAM: 16 GB, 2 x 8 GB, Disk: 1 TB SSD GPU : RTX 3060 6 GB		✓	✓
Legion Pro 5i Gen 8 Intel (16") CPU: i7-13700HX RAM: 16 GB Disk: 512 GB SSD GPU: RT 4060 8GB		✓	✓
Legion Pro 7i Gen 8 Intel (16") with RTX 4070 CPU: i9-13900HX RAM:16GB Disk: 512 GB SSD GPU : RTX 4070 8GB		✓	✓

For more details, please visit:
https://www.dell.com/en-us/shop/dell-laptops/sr/laptops （Dell）
https://www.lenovo.com/us/en/d/gaming-laptops/?IPromoID=LEN344641（Legion）

Panda Scanner

www.panda-scanner.com info@panda-scanner.com +86 512 6731 7873